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INSTRUCTIONS FOR USE - NEW OSTEO

NEW OSTEO - CALCIUM SULPHATE BONE VOID FILLER

1) PRODUCT DESCRIPTION AND MATERIALS

NEW OSTEO Calcium Sulphate Bone Void Filler pellets and powder are provided sterile for single patient use. The biodegradable, radiopaque pellets are resorbed in approximately 60 days when used in accordance with the device labelling. Numerous bone graft operations are undertaken every year to repair bone defects caused by surgery, tumours, trauma, infections and other damage. NEW OSTEO Calcium Sulphate Bone Void Filler is intended for use in clinical situations where the use of autologous grafts or other bone graft substitutes may be undesirable, due either to the risk of associated infection or unavailability. NEW OSTEO is manufactured from medical grade calcium sulphate dihydrate ($\text{CaSO}_4 \cdot 2\text{H}_2\text{O}$) of high purity, that resorbs and is replaced with bone during the healing process. Also, as the bone void filler pellets and powder are biodegradable and biocompatible, they may be used at an infected site associated with antibiotic as delivery device. Either the pellets or the powder can be mixed with autologous growth factor from platelets concentration supplied by SEFC machine (Biológica - GMReis Company) resulting in an osteoconductive and osteoinductive gel to enhance bone healing.

2) INTENDED USES

NEW OSTEO is indicated only for bone voids or defects / gaps that are not intrinsic to the stability of the bone structure. NEW OSTEO is indicated to be gently packed into bone voids or defects of the skeletal system (i.e., long bones, extremities, spine and pelvis). Pellets should not be crushed shaved or cut and powder must be mixed to the saline solution to hydrate, molded and packed into non-bearing voids and gaps. The bone defects or cavities may be surgically created or as a result of traumatic injury. NEW OSTEO provides a bone graft substitute that resorbs and is replaced with bone during the healing process. NEW OSTEO is biodegradable, biocompatible and radiopaque and may be used at any infected site.

3) INDICATIONS FOR USE

NEW OSTEO is indicated to fill a bone void or defect created by:

- Surgery;
- Cyst;
- Tumors;
- Osteomyelitis or other traumatic injury;

4) CONTRAINDICATIONS

NEW OSTEO isn't indicated for:

- Filling of defects which are intrinsic to the stability of the bone structure;*
- Defects which require the pellets to be crushed;
- Severe vascular or neurological disease;
- Uncontrolled diabetes;
- Severe degenerative bone disease;
- Pregnancy;
- Uncooperative patient who can't or won't follow postoperative instructions including individuals who abuse drugs or alcohol;
- Hypercalcaemia;

* NEW OSTEO can not be used to support mechanical loads.

5) INSTRUCTIONS FOR USE

1. Accurate diagnosis of any bone cyst in an essential pre-requisite
2. Expose the bone defect. Raise a window of cortical bone with attached muscle or soft tissue when possible. Alternatively raise periosteum without dissection in the extraperiosteal plane. Decortication is to be preferred by sharp osteotome to raise small slivers of bone attached to the deep surface of periosteum.
3. In case of a cyst, remove with a curette or preferably a high-speed burr. A dental mirror or arthroscope improves clearance of all unwanted tissue by improving visualization.
4. Lavage the prepared bone graft site with a saline solution until a healthy bleeding surface is developed on all sides. Pellets and molded implant from powder should not be crushed, shaved or cut.
5. Determine the size of the bone void by filling with the saline solution syringe. The bone defect should then be tightly packed with NEW OSTEO.
6. Cover the opening of the bone void with the preserved cortical bone window or flap. Alternatively suture the periosteal flap firmly across the pellets or molded implant of NEW OSTEO.

It is of benefit to create a separate micro environment deep to the periosteum. Clean the wound in layers without tension.

6) PRECAUTIONS

The patient should be advised to report any related pain, decrease in range of motion, swelling, fever and unusual incidences. The patient is to be cautioned to govern activities and protect the surgical site from unreasonable stresses, and follow the instructions of the physician with respect to follow-up care and treatment. The patient is to be warned of surgical risks, and made aware of possible adverse effects.

7) ADVERSE EFFECTS

Peripheral neuropathies have been reported following surgery. Subclinical nerve damage occurs more frequently, possibly the results of surgical trauma. Material sensitivity/allergic reactions in patients following surgery have been rarely reported. Implantation of foreign material in tissues can result in histological reactions involving macrophages and fibroblasts. The clinical significance of this effects is uncertain, as similar changes may occur with the precursor or during the healing process. Infection can lead to failure / removal of the void filler.

8) STERILIZATION

Unless opened or damaged, pellets and powder and mixing solution are supplied sterile in jars that are double pouched. Check all packaging for punctures or other damage. Sterilization of components is by gamma radiation. When removing the NEW OSTEO, appropriate aseptic procedure must be observed. Non-aseptic personnel may open the outer pouch. Aseptic personnel must retrieve the inner barrier and the pellet jar. Resterilisation of the NEW OSTEO by steam autoclave or ethylene oxide gas is prohibited, because it will lose its properties.

9) STORAGE CONDITIONS

NEW OSTEO must be stored and kept in a clean room at temperatures between 4 and 45°C.

10) PRODUCT COMPLAINTS

Any Health Care Professionals, e.g., customer users of GMReis products, who have complaints or have experienced dissatisfaction in the product, quality, identify, durability, reliability, safety, effectiveness, and/or performance should notify the GMReis distributor or GMReis directly by e-mail qualidade@gmreis.com.br or by below address. When filing a complaint, please note the component name and any other descriptive information, plus your name, address, and nature of complaint.

GM dos Reis Ind. e Com. Ltda.
600, Lote 3 Quadra F - Pierre Simon de Laplace Ave
Techno Park
Campinas / SP - Zip Code 13069-320 - Brazil
Phone: + 55 19 3765-9900 Fax: + 55 19 3765-9111
E-mail: qualidade@gmreis.com.br
www.gmreis.com.br